

CardioTwin AI — Product Requirements Document

Version: 1.0 | **Date:** February 21, 2026 | **Status:** FINAL — Ready for 24-Hour Sprint

1. Mission & Theme Alignment

Mission: Build a working prototype in 24 hours that demonstrates real-time cardiometabolic risk scoring from live sensor data, delivered through a web dashboard and WhatsApp notifications.

Hackathon Theme: "From Data to Prevention: AI as Your Health Partner"

Theme Word	Our Implementation
Data	Live biometric data from MAX30102 + NTC sensors
Prevention	Predictive risk scoring + what-if projections
AI	Adaptive risk engine + pattern detection
Health Partner	Real-time nudges via WhatsApp + dashboard

2. Problem Statement

- **~11% of all deaths** in Nigeria are attributed to CVDs (WHO)
- **70% of cardiovascular deaths** occur in people who were **never screened**
- Nigeria has **~4 physicians per 10,000 people** — preventive cardiology is nearly non-existent

- No existing system makes the connection between "**what I just did**" and "**what my heart just felt**" immediate, tangible, and personal

Our Hypothesis: If we give people a real-time, personalized score that visibly reacts to their physiological state, they will make better health decisions because consequences become immediate — not abstract.

3. Solution Overview

CardioTwin AI is a station-based health screening system with three layers:

1. **Sense** — Captures HR, HRV, SpO₂, and skin temperature from a finger-placed sensor station
2. **Score** — AI risk engine processes raw biometrics into a **CardioTwin Score (0-100)**
3. **Nudge** — Delivers alerts via web dashboard, WhatsApp, and physical vibration buzz

Scope: IN □

- Working hardware station with finger sensor
- Real-time data streaming to Azure cloud
- CardioTwin Score algorithm (0-100, 4 zones)
- Live web dashboard with graphs and score gauge
- WhatsApp alert integration via Twilio
- What-if risk projection (static/pre-calculated)
- Before/after demo capability
- Vibration motor alert on station

Scope: OUT □

- Wearable device (pitched as v2 vision only)
- Facial expression analysis

- Mobile native app (web PWA only)
 - User authentication / multi-user accounts
 - Deep learning models
 - Medical-grade diagnostic claims
-

4. System Architecture

CARDIOTWIN AI

LAYER 1: SENSING

HEALTH STATION

ESP32
+ MAX30102
+ NTC Therm.
+ Vibration
Motor

WiFi / HTTP POST


LAYER 2: MODELING

AZURE CLOUD

Data Ingestion
(FastAPI)

Adaptive Risk
Engine
(Scikit-learn)

CardioTwin
Score API
+ Session DB



5. The CardioTwin Score

The single most important feature — transforms complex multi-dimensional biometric data into one number anyone can understand.

Score Zones

Zone	Range	Label	Trigger
🟢 Green	80-100	Thriving	No alert
🟡 Yellow	55-79	Mild Strain	WhatsApp nudge
🟠 Orange	30-54	Elevated Risk	WhatsApp + vibration buzz
🔴 Red	0-29	Critical Strain	WhatsApp + vibration + urgent message

Scoring Weights

Component	Weight	Why
HRV (RMSSD)	40%	Strongest non-invasive predictor of cardiac autonomic health
Heart Rate	25%	Direct indicator of cardiovascular load
SpO ₂	20%	Blood oxygenation — safety check for respiratory distress
Skin Temperature	15%	Correlates with stress response and inflammation

How It Works

1. **Baseline Calibration (first 30s):** Records resting HR, HRV, SpO₂, and temperature as personal reference
2. **Real-Time Scoring (every 2s):** Compares incoming data against baseline, calculates weighted composite
3. **Pattern Detection (over sessions):** Identifies personal risk patterns over time

6. Hardware Specification

Station Design (Not a wearable — by design)




The MAX30102 requires firm, stable skin contact and minimal motion for accurate PPG readings. A station-based finger sensor gives **clinical-grade accuracy** every time, eliminates motion artifact risk, and avoids demo failures.

During the pitch, frame the wearable as the v2 vision: "The algorithm is built. The form factor evolves."

Components

Component	Purpose	Est. Cost
ESP32-WROOM-32	Microcontroller + WiFi	~ ₦ 3,500
MAX30102	Heart Rate, HRV, SpO ₂	~ ₦ 2,000
10kΩ NTC Thermistor	Skin Temperature	~ ₦ 300
10kΩ Resistor	Voltage divider for NTC	~ ₦ 50
Mini Vibration Motor (3V)	Tactile alert	~ ₦ 500
2N2222 NPN Transistor + 1kΩ Resistor	Motor driver	~ ₦ 100
Breadboard + Jumper Wires	Assembly	~ ₦ 800
Micro-USB Cable	Power + programming	~ ₦ 500
TOTAL		~₦7,750 (~\$10)

Wiring Reference

ESP32 PIN	→	COMPONENT
3.3V	→	MAX30102 VIN
GND	→	MAX30102 GND
GPIO 21 (SDA)	→	MAX30102 SDA
GPIO 22 (SCL)	→	MAX30102 SCL
3.3V	→	10kΩ Resistor 
GPIO 34 (ADC)	→	Junction 
GND	→	NTC Thermistor 
GPIO 25	→	1kΩ → 2N2222 Base Collector → Motor → 3.3V Emitter → GND

Station Enclosure

Simple cardboard/foam box with a hole for the MAX30102 sensor on top. Labeled "Place Finger Here 🖐️". Sensor MUST be secured firmly (hot glue) to prevent wobble during demo.

7. API Contract

All endpoints served from: `https://cardiotwin.azurewebsites.net`

POST /api/session/start

Starts a new measurement session.

Request:

```
{ "session_id": "demo", "user_phone": "+2348012345678" }
```

Response:

```
{ "status": "session_started", "session_id": "demo" }
```

POST /api/reading

Receives biometric reading from ESP32. Called every 2 seconds.

Request:

```
{  
  "bpm": 72, "hrv": 42.3, "spo2": 98.1,  
  "temperature": 36.4, "timestamp": 45000, "session_id": "demo"  
}
```

Response (calibrating):

```
{  
  "status": "calibrating", "readings_collected": 8,  
  "readings_needed": 15, "alert": false  
}
```

Response (scored):

```
{  
  "status": "scored", "score": 86.2, "zone": "GREEN",  
  "zone_label": "Thriving", "zone_emoji": "🌱",  
  "alert": false, "nudge_sent": false,  
  "components": {  
    "heart_rate": { "value": 72, "score": 95.2 },  
    "hrv": { "value": 42.3, "score": 88.1 },  
    "spo2": { "value": 98.1, "score": 100.0 },  
    "temperature": { "value": 36.4, "score": 93.5 }  
  },  
  "baseline": {  
    "resting_bpm": 71.5, "resting_hrv": 43.1,  
    "normal_spo2": 98.0, "normal_temp": 36.35  
  }  
}
```

GET /api/score/{session_id}

Returns latest score for frontend polling.

GET /api/history/{session_id}

Returns all scores for chart rendering. Array of score objects with timestamps.

POST /api/predict

What-if risk projection.

Request:

```
{ "session_id": "demo", "days": 90 }
```


Response:

```
{  
  "current_score": 41.0, "projected_score": 35.2,  
  "projected_resting_hr_increase_bpm": 6.8,  
  "current_risk_category": "Elevated Risk",  
  "projected_risk_category": "Critical Strain",  
  "disclaimer": "Statistical projection only. Not a medical diagnosis."  
}
```

8. Tech Stack

Layer	Technology	Source
Microcontroller	ESP32	Off-the-shelf
Sensors	MAX30102 + NTC Thermistor	Off-the-shelf
Firmware	Arduino C++	Open source
Backend	Python FastAPI	Open source
AI/ML	Scikit-learn, NumPy	Open source
Cloud	Azure App Service + Cosmos DB	GitHub Student Dev Pack
Frontend	Next.js 14 + Tailwind CSS + Recharts	Open source
Notifications	Twilio WhatsApp + SMS	GitHub Student Dev Pack
Version Control	GitHub + GitHub Copilot	GitHub Student Dev Pack

9. Team Roles & Responsibilities



Person 1 — Hardware Engineer

Deliverable	Priority
Wire ESP32 + MAX30102 + NTC + vibration motor	☐ Critical
Write firmware: read HR, HRV, SpO ₂ , temperature	☐ Critical
Implement WiFi + HTTP POST to backend every 2s	☐ Critical
Vibration motor trigger on API alert response	☐ Medium
Clean station enclosure with finger label	☐ Medium
Calibrate sensor accuracy + integration test	☐ Critical

Key Specs: - Send JSON via HTTP POST every 2 seconds - Finger detection: IR > 50,000 = finger present - Reject HR outside 30-220 BPM - HRV = RMSSD from beat-to-beat intervals - SpO₂ from RED/IR ratio (linear approximation) - Vibrate when API response contains `"alert": true`

Person 2 — AI Engineer

| Deliverable | Priority | |---|---|---| | Design CardioTwin Score algorithm (weighted composite) | ☐ Critical | | Implement baseline calibration (first 15 readings) | ☐ Critical | | Build zone classification (Green/Yellow/Orange/Red) | ☐ Critical | | Build anomaly detection for real-time alerts | ☐ Medium | | Build what-if projection engine (trend extrapolation) | ☐ Medium | | Write nudge message templates per zone | ☐ Critical | | Tune weights/thresholds during integration | ☐ Medium |

Person 3 — Backend Developer (Python)

Deliverable	Priority
Scaffold FastAPI project, deploy to Azure	<input type="checkbox"/> Critical
Implement <code>/api/reading</code> endpoint	<input type="checkbox"/> Critical
Implement <code>/api/score</code> + <code>/api/history</code> endpoints	<input type="checkbox"/> Critical
Implement <code>/api/session/start</code> endpoint	<input type="checkbox"/> Critical
Integrate Person 2's scoring engine	<input type="checkbox"/> Critical
Set up Twilio WhatsApp integration	<input type="checkbox"/> Critical
Implement <code>/api/predict</code> endpoint	<input type="checkbox"/> Medium
Set up CORS for frontend access	<input type="checkbox"/> Critical

Person 4 — Frontend Developer

Deliverable	Priority
Scaffold Next.js + Tailwind CSS project	<input type="checkbox"/> Critical
Build animated CardioTwin Score gauge (circular, color-coded)	<input type="checkbox"/> Critical
Build real-time biometric cards (HR, HRV, SpO ₂ , Temp)	<input type="checkbox"/> Critical
Build real-time score line chart (Recharts)	<input type="checkbox"/> Critical
Implement 2-second polling from backend API	<input type="checkbox"/> Critical
Build before/after comparison view	<input type="checkbox"/> Medium
Build what-if simulator panel	<input type="checkbox"/> Medium
Polish: animations, transitions, responsive design	<input type="checkbox"/> Medium

Person 5 — Pitch Lead & PM

Deliverable	Priority
Write pitch script (3 minutes)	<input type="checkbox"/> Critical
Build slide deck (8 slides max)	<input type="checkbox"/> Critical
Coordinate integration between all members	<input type="checkbox"/> Critical
Record backup demo video	<input type="checkbox"/> Critical
Run 3+ dress rehearsals	<input type="checkbox"/> Critical
Prepare Q&A answers for judges	<input type="checkbox"/> Critical
Research Nigeria CVD stats for pitch	<input type="checkbox"/> Medium

10. 24-Hour Sprint Timeline

HOUR	TASK	OWNER
0-1	Kickoff: scope lock, repo setup, Azure provisioning, agree on API contracts	ALL
1-4	Hardware: wire + test raw sensor readings AI: design + code score algorithm Backend: scaffold FastAPI + deploy Azure Frontend: scaffold Next.js + design UI PM: draft pitch narrative + start deck	Person 1 Person 2 Person 3 Person 4 Person 5
4-8	Hardware: implement WiFi streaming AI: baseline calibration + zone logic Backend: /api/reading + /api/score live Frontend: score gauge + biometric cards PM: coordinate, write demo script	Person 1 Person 2 Person 3 Person 4 Person 5
8-12	Hardware: test continuous streaming AI: what-if projection engine Backend: Twilio WhatsApp + /api/predict Frontend: live chart + what-if panel PM: first end-to-end integration test	Person 1 Person 2 Person 3 Person 4 Person 5
12-14	☐ MANDATORY REST – 2 HOURS (You WILL make worse decisions without rest)	ALL
14-18	Full integration: hardware→backend→frontend Bug fixing + edge case handling PM: record backup demo video	ALL ALL Person 5
18-21	Polish: UI animations, loading states Calibrate score thresholds with real data Build clean station enclosure Finalize pitch deck	Person 4 Person 2 Person 1 Person 5
21-23	Full dress rehearsal (3 TIMES MINIMUM) Fix any demo-breaking bugs	ALL ALL

23-24 Final prep: charge devices, test venue
WiFi, load backup video, deep breaths

ALL

11. Demo Script (3 Minutes)

Act 1 — The Hook (30s)

"Every 33 seconds, someone dies of cardiovascular disease. In Nigeria, 7 out of 10 of those people were never screened. What if your body could warn you — not your doctor, not a hospital — but your body, in real time, before it's too late?"

Act 2 — Live Demo (90s)

1. Volunteer places finger on CardioTwin Station
2. Dashboard shows data streaming live → Score: **86 — Thriving**
3. Volunteer does 30s of burpees (presenter narrates the AI engine)
4. Volunteer returns, places finger again
5. Score drops live: **41 — Elevated Risk**
6. Vibration motor buzzes on station
7. WhatsApp message arrives on volunteer's phone
8. Show What-If: "If this pattern continues 90 days → risk category shifts"

Act 3 — The Vision (45s)

"Today: one person, one station, one score. Tomorrow: community health workers across Nigeria, each with a ₦8,000 CardioTwin device, screening entire villages in hours. Population-level dashboards helping state health ministries allocate resources where risk is highest."

Act 4 — The Close (15s)

"From data, to prevention, to partnership. This is CardioTwin AI — your heart's early warning system."

□ BACKUP PLAN

Pre-record a full demo video the night before. If hardware fails, play the video. **Never demo without a backup.**

12. Risk Register

Risk	Likelihood	Impact	Mitigation
WiFi fails at venue	High	Critical	Test venue WiFi early. Backup: mobile hotspot. Emergency: pre-recorded video
MAX30102 gives erratic readings	Medium	High	Station-based (not wearable) = stable contact. Calibrate before demo. Have backup recorded data
Azure deployment issues	Medium	High	Deploy by Hour 8. Test early. Fallback: run locally on laptop
Twilio WhatsApp fails	Medium	Medium	Pre-join Twilio sandbox. Test with real number. Fallback: show SMS or mock notification
Team member burns out	High	High	Mandatory 2-hour rest at Hour 12-14. No hero coding
Integration breaks at Hour 20	Medium	Critical	Integrate early (Hour 8 first test). API contracts agreed at Hour 0
Judges question "Digital Twin" claim	High	High	We DON'T call it a digital twin. It's an "Adaptive Risk Model" — honest and defensible

13. Prepared Judge Q&A

Q: "Is this medically accurate?"

"CardioTwin is a wellness screening tool, not a medical device. Our MAX30102 sensor provides research-grade PPG data, and our scoring algorithm is based on published clinical correlations between HRV and cardiovascular risk. We clearly disclaim this in the UI."

Q: "Why not a wearable?"

"Clinical-grade PPG requires stable skin contact. Even Apple Watch uses motion compensation algorithms that took years to develop. Our station gives accurate data every time. The wearable is our v2 — the intelligence is built today."

Q: "How is this different from a Fitbit?"

"Fitbit gives you data. We give you a decision. Our CardioTwin Score collapses complex biometrics into one actionable number, and our nudge system delivers guidance through WhatsApp — where 93 million Nigerians already live. No app download required."

Q: "Can this scale?"

"The hardware costs ₦8,000. The backend runs on Azure free tier. Community health workers could screen a village of 500 in a single day. The data aggregates into population-level risk heatmaps for state health ministries."

14. Data Privacy & Ethics

Principle	Implementation
Encryption	TLS 1.3 in transit, Azure managed keys at rest
Compliance	Follows NDPR (Nigeria Data Protection Regulation)
Consent	Explicit opt-in before data collection
No Diagnosis	Wellness tool — clearly disclaimed throughout
Data Ownership	Users can export or delete all data

15. Business Model (For Pitch)

Stream	Model
B2C	Freemium — free basic score, paid insights (₦2,000/mo)
B2B	Corporate wellness dashboards, HMO partnerships (Leadway, AXA Mansard)
B2G	State health ministry screening tools, population CVD heatmaps

16. What Makes This Win

Differentiator	Why It Matters
Working hardware prototype	90% of teams are software-only
CardioTwin Score (0-100)	One number anyone understands
WhatsApp delivery	Shows deep understanding of Nigerian users
\$10 hardware cost	Proves accessibility at scale
Honest AI claims	"Adaptive Risk Model" builds trust with judges
Clear business model	Shows maturity beyond a hackathon
Backup demo video	Professional risk management

Remember: You're not building a product in 24 hours. You're building a story backed by a working prototype. The story wins. The prototype proves you can deliver.

****Go build it. ☐☐☐**